


AlumaSeal CS Sealing films

Specially formulated for cold storage to -80 °C
Excellent barrier properties to retard evaporation
Easily pierceable with pipet tips or robotic probes for sample recovery
Not recommended for PCR or thermalcycling

Specification :

- ✓ Material : aluminum
- ✓ Foil thickness : 50 µm
- ✓ Adhesive material : acrylic
- ✓ Non-perforated end tab (9.5mm)
- ✓ Dimensions : 82.6 x 132.6 mm
- ✓ Temperature range -80°C to +120°C
- ✓ DNase-, RNase-, and nucleic-acid-free
- ✓ Sterile



Unités :
50 films