

<b>Specification#: QA-FFGEL0030</b>	<b>Date of Revision:</b> 08/09/16	<b>Revision:</b> 6
	<b>Original Date:</b> 13/02/12	
<b>Item Number: FFGEL0030</b>	<b>Reason for Revision:</b> Tolerance and palletisation modifications	
	<b>Originator:</b> Tixier Angeline	
	<b>Approved:</b> Boulay William	<b>Date:</b> 08/09/16

## General Product Information

Product Name	Nylon Gel 250 g
Description	Nylon Gel 250g are composed of a sealed Polyethylene sack that is filled with eutectic gel. It is used for keeping products cold inside an isothermal packaging.
Item Number	FFGEL0030

## Finished Product Components

Pouch component	Polyethylene
Cryopak Solution	Gels made with potable Water





## Specifications

Outer Dimensions (with weldings)	160 x 155 x 20 mm $\pm 10\%$
Weight	250g $\pm 5\%$
Units per carton	36
Units per pallet	1620
Visual Appearance	See Diagram 1

## Physical Properties

Freezing point	32°F (0°C)
MSDS	Available on Request

## Usage Warning

-  Do not stack pallets ;
-  Do not puncture or tear the gel pack (discard if gel pack leaks) ;
-  Open carefully the carton box (without a craft knife) ;
-  Always freeze the gel pack in flat position.

## Company Information

Manufacturer's name : TCP France  
 Street Address : 1 rue Berthelot  
 City, State, Zip code : Maromme, 76150 France  
 Emergency telephone no. : +33 (0) 232825965

# Specification Sheet - External

<b>Specification#:</b> QA-FFGEL0030	<b>Date of Revision:</b> 08/09/16	<b>Revision:</b> 6
	<b>Original Date:</b> 13/02/12	
<b>Item Number:</b> FFGEL0030	<b>Reason for Revision:</b> Tolerance and palletisation modifications	
	<b>Originator:</b> Tixier Angeline	
	<b>Approved:</b> Boulay William	<b>Date:</b> 08/09/16

Diagram 1 (may not be actual size)

