



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Germany
+49 (0)6221 / 8442-44
<https://www.neolab.de/en>

VAT No.: DE 143 450 657



neoBox-200, storage box | 2-1968

Product Images



Description

Stackable storage box made from PP with 20 x 10 compartments for vessel up to 12 mm Ø, maximum height 46 mm (HPLC vials, cryo tubes, reaction vessels 1.5/2.0 ml).

Tight closing lid with hinges.

Alphanumeric grid on the lid.

Temperature-resistant between -80°C and 121°C

Material: polypropylene, autoclavable up to 121°C

Dimensions: 275 x 140 x 60 mm (W x D x H)

Accessory: Insert for reaction vessels 1.5 ml (Order No. 2-1210)

DOMINIQUE DUTSCHER SAS

Additional Information

For small bottle/ vial type	vials, 12 mm Ø
Storage type	boxes
Alphanumeric coding	Yes
Material	polypropylene (PP)
Closure/ lid type	Klappverschluss/-deckel
With special chemical resilience	No
Height	60 mm
Length	140 mm
Width / Breadth	275 mm
Color	transparent (milky)
Stackable	Yes
Suited for	tubes/ vials microvessels/ vials
Autoclavable	Yes
With closure/ lid	Yes
Suited for cooling in	Ultratiefkühlschrank/ Deep Freezer
For containers with diameter up to	12 mm
Capacity	200
For containers with volume up to	1,2 ml 2 ml
For small tube/ micro vessel type	reaction vessels cryogenic tubes
With labeling field	No
Closure/ lid available as accessory	No
In a set	No
With beta-radiation protection	No
For different vessel types/ volumes	Yes
Wells (W x L)	20 x 10
For light sensitive samples	No
With removable grid inlet	Yes

With grit divider	Yes
Grid dividers available as accessories	Yes
Available in special sizes	No
Shipping conditions	Raumtemperatur

DOMINIQUE DUTSCHER SAS