

TECHNICAL DATA SHEET

Part number	764047
Description	Tube rack stackable



Description :

- ✓ Storage type : racks
- ✓ Alphanumeric coding : Yes
- ✓ Material : polypropylene (PP)
- ✓ With special chemical resilience : Yes
- ✓ Height : 70 mm
- ✓ Length : 146 mm
- ✓ Width / Breadth : 146 mm
- ✓ Color : orange
- ✓ Stackable : Yes
- ✓ Suited for : variable applications
- ✓ Autoclavable : Yes
- ✓ Suited for cooling in: freezer
- ✓ For containers with diameter up to : 13 mm
- ✓ Capacity : 81
- ✓ For small tube/ micro vessel type : test tubes|centrifuge tubes
- ✓ Wells (W x L) : 9 x 9
- ✓ Compartment size : 13 x 13 mm
- ✓ Suited for incubators Yes

Pieces by sales :

1 Rack
