

ClearLine AeraSeal sealing films

Application :


Cell, tissue, bacterial culture

Advantage :


Gas permeable

Compatible with :

Skirted and semi-skirted


Features

Packaging	Roll
End-tabs	Features two end-tabs
Collar	Without collar
Material	Rayon
Adhesive material	Acrylic
Piercable	Yes
Chemically resistant	No
Breathable	Yes
Opacity	-
Colour	-
Sterile	Yes
Functional temperature range	-20°C to +80°C
Film/Foil thickness	139,7 µm
Adhesive thickness	-
Dimensions (without end-tabs)	120,7 x 82,6


Units : 2 rolls of 50 films