

### ClearLine AeraSeal sealing films

**Application :**


Cell, tissue, bacterial culture

**Advantage :**


Gas permeable

**Compatible with :**

Skirted and semi-skirted


**Features**

<b>Packaging</b>	Roll
<b>End-tabs</b>	Features two end-tabs
<b>Collar</b>	Without collar
<b>Material</b>	Rayon
<b>Adhesive material</b>	Acrylic
<b>Piercable</b>	Yes
<b>Chemically resistant</b>	No
<b>Breathable</b>	Yes
<b>Opacity</b>	-
<b>Colour</b>	-
<b>Sterile</b>	No
<b>Functional temperature range</b>	-20°C to +80°C
<b>Film/Foil thickness</b>	139,7 µm
<b>Adhesive thickness</b>	-
<b>Dimensions (without end-tabs)</b>	120,7 x 82,6



**Units : Starter kit with SealMate dispenser and 2 rolls of 50 films**