

ClearLine AeraSeal sealing films

Application :


Cell, tissue, bacterial culture

Advantage :


Gas permeable

Compatible with :

Skirted and semi-skirted


Features

| | |
|--------------------------------------|-----------------------|
| Packaging | Roll |
| End-tabs | Features two end-tabs |
| Collar | Without collar |
| Material | Rayon |
| Adhesive material | Acrylic |
| Piercable | Yes |
| Chemically resistant | No |
| Breathable | Yes |
| Opacity | - |
| Colour | - |
| Sterile | Yes |
| Functional temperature range | -20°C to +80°C |
| Film/Foil thickness | 139,7 µm |
| Adhesive thickness | - |
| Dimensions (without end-tabs) | 120,7 x 82,6 |



Units : Starter kit with SealMate dispenser and 2 rolls of 50 films