

ClearLine ThermalSeal RTS sealing films

Application :


qPCR, Long-term storage, Crystallography, HPLC

Advantage :


The encapsulated silicone adhesive is non-tacky until pressed against the sealing surface

Compatible with :

All kind of skirts


Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	With or without collar
Material	Polyolefin
Adhesive material	Silicone
Piercable	No
Chemically resistant	Yes
Breathable	No
Opacity	Optically Clear
Colour	-
Sterile	Yes - DNase RNase free
Functional temperature range	-80°C to +105°C
Film/Foil thickness	50,8 µm
Adhesive thickness	50,8 µm
Dimensions (without end-tabs)	116,2 x 78,3


Units : 50 films