

## ClearLine Ezcap 1x8 PCR FoilSeal sealing films

**Application :**


PCR

**Advantage :**


Contains 6 strips on a continuous liner

**Compatible with :**

8-Tube Strip and Strip wells - Skirted and semi-skirted


**Features**

<b>Packaging</b>	Films
<b>End-tabs</b>	Features two end-tabs with perforations
<b>Collar</b>	Without collar
<b>Material</b>	Aluminium
<b>Adhesive material</b>	Acrylic
<b>Piercable</b>	Yes
<b>Chemically resistant</b>	No
<b>Breathable</b>	No
<b>Opacity</b>	Light blocking
<b>Colour</b>	-
<b>Sterile</b>	No - DNase RNase free
<b>Functional temperature range</b>	-20°C to +120°C
<b>Film/Foil thickness</b>	50,8 µm
<b>Adhesive thickness</b>	36,8 µm
<b>Dimensions (without end-tabs)</b>	81,2 x 10,1


**Units : 300 films**