



Ampules, m-ColiBlue24 broth, glass , pk/20



Réf. produit: 2608420

Read and confirm results in 24 hours

Get superior sensitivity - 1 CFU/100mL

Conduct emergency testing and routine monitoring

Eliminate special equipment and ultraviolet lamp

Minimize background growth of non-coliform bacteria

▼ Caractéristiques

Matériel:	Glass
Media Type:	m-ColiBlue24 Broth
Shelf Life (days) as of shipment date:	150
Volume / taille du carton:	20 /pk