

Vegetable LB Broth (Lennox)

Cat. 2038

Recommended medium for maintaining and cultivating recombinant strains of E. coli.

Practical information

Applications	Categories
Preparation and recovery of competent cells	Escherichia coli
Industry: Culture media for Molecular biology	



Principles and uses

Vegetable LB Broth (Lennox) is a nutritionally rich medium based on the LB Broth (Lennox) (Cat. 1163) designed as an alternative to classical animal-based media for the growth and maintenance of pure cultures of recombinant strains of E. coli used in molecular microbiology procedures.

These strains are generally derived from E. coli K12, which are unable to produce vitamin B, so this media is formulated to enhance the growth of nutritionally demanding microorganisms. This strain of E. coli has been further modified through specific mutation to create an auxotrophic strain that is not capable of growth on nutritionally deficient media.

In this case, the ingredients are animal-free in order to minimize the risk of bovine spongiform encephalopathy in culture media containing bovine materials. Vegetable origin peptone provides nitrogen, vitamins, minerals and amino acids essential for growth. Yeast extract is source of vitamins, particularly the B-group essential for bacterial growth. Sodium chloride supplies essential electrolytes for transport and osmotic balance.

Vegetable LB Broth (Lennox) has a different sodium chloride level than other media such as Luria Broth (Miller LB Broth) (Cat. 1551) or Luria Broth (Miller's Modification) (Cat. 1266). This allows to select the optimum salt concentration of the medium for a specific strain.

Formula in g/L

Sodium chloride	5	Yeast extract	5
Vegetable peptone	10		

Preparation

Suspend 20 grams of the medium in one liter of distilled water. Mix well and dissolve by heating with frequent agitation. Boil for one minute until complete dissolution and dispense into tubes. Sterilize in autoclave at 121 °C for 15 minutes.

Instructions for use

- Carry out the experimental procedure according to appropriate use or purpose.
- Inoculate and incubate at a temperature of 35±2 °C for 18-24 hours.

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
w/o rests	Fine powder	Beige	Clear beige	7,0±0,2

Microbiological test

Incubation conditions: (35±2 °C / 18-24 h).

Microorganisms	Specification
Escherichia coli ATCC 23724	Good growth
Escherichia coli ATCC 33694	Good growth
Escherichia coli ATCC 33849	Good growth
Escherichia coli ATCC 39403	Good growth
Escherichia coli ATCC 47014	Good growth

Storage

Temp. Min.:2 °C
Temp. Max.:25 °C

Bibliography

- Ausubel, Brent, Kingston, Moore, Seidman, Smith and Struhl (ed.). 1994. Current protocols in molecular biology, vol. 1. Green Publishing Associates, Inc., Brooklyn, N.Y.
- Mackenzie E.F.N., Taylor E.N. and Gilbert W.E. 1948. Recent experiences in the rapid identification of E. coli, type I. J. Gen. Microbiol. 2:197-204.
- Atlas, R.M., L.C. Parks (1993) Handbook of Microbiological Media. CRC Press, Inc. London
- Lennox. (1955). Virology 1:190.
- Sambrook, Fritsch and Maniatis. (1989). Molecular cloning: a laboratory manual, 2nd ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y.