

Bacteriological Ox Bile

Cat. 1710

Ingredients (Biles)

Practical information

Principles and uses

Bacteriological Ox Bile is prepared by a low temperature dehydration process. It is used as a selective inhibitory agent in culture media such as Brilliant Green Bile Broth 2%.

Physical-chemical characteristics

Description	Specification	Typical Analysis
Loss on drying	<6%	3,1%
pH (2% solution)	7,5-9,2	8,3
Colic acid	>45%	47%

Growth supporting properties

Descripción	Value
Brilliant Green 2% Bile Broth	Good/Bueno

Microbiological test

Description	Specification
Hongos y levaduras	100 CFU/g
Coliformes	Negative
Salmonella	Negative
Recuento en placa	<5.000 CFU/g

Storage

Temp. Min.:2 °C
Temp. Max.:25 °C