



PRACTICAL INNOVATION

Ziath

Tube Picker

Simplifies retrieving cryogenically frozen samples

Ziath's Tube Picker allows one-handed removal of microtubes from racks.



The Tube Picker simplifies retrieval of tubes, minimises heat transfer and reduces the risk of cross contamination.



Features

- Single handed operation
- Tube selection from anywhere in the rack
- For use with 8.5mm tubes
- Easy-grip
- Reduces risk of thawing and contamination
- Simplifies retrieving samples
- Compatible with Matrix, Nunc, FluidX, ABgene, REM-P, Micronic, and most other microtubes

Ordering Information

Code	Description	Quantity
ZPIK-01	Tube Picker	1