

Data sheet

12.10.2020

Semi-quantitative test strips QUANTOFIX Phosphate

Parameter*	Value*
REF	91320
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Phosphate
Gradation	0 • 3 • 10 • 25 • 50 • 100 mg/L PO ₄ ³⁻
Method	Phosphomolybdenum blue
Color change	white → blue-green
Instrumental measuring range	3–80 mg/L PO ₄ ³⁻

*These data are typical values and must not be regarded as specification.

www.mn-net.com



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE/international:
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
E-mail: info@mn-net.com

CH:
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
E-mail: sales-ch@mn-net.com

FR:
Tel.: +33 388 68 22 68
Fax: +33 388 51 76 88
E-mail: sales-fr@mn-net.com

US:
Tel.: +1 484 821 0984
Fax: +1 484 821 1272
E-mail: sales-us@mn-net.com