

## Description

- SI Analytics piston burette TITRONIC® 500 with magnetic stirrer and 20 ml exchangeable head

## Part number

- 780443

Pieces by sales : 1



---

## Description

- Motor driven piston burette for exchangeable heads for dosing and titration
- The TITRONIC® 500 is the ideal piston burette for manual titrations, accurate dosing applications as well as the preparation of solutions. When used with TitrSoft 3.3 and TitroLine® 7000, 7750 and 7800 it acts as an automatic dispenser or a titration burette.
- General features
  - High visibility, full color display that can be easily viewed from a distance and at extreme angles.
  - Reagent data is securely stored in the intelligent and interchangeable modules
  - Includes up to three USB, one LAN and two RS232 ports for expansion and connection of devices
  - Export the results as PDF or CSV, also to networks.
  - User management/password
- Specific features for the TITRONIC ® 500
  - Manual titration with the hand controller (mouse)
  - Titration rate can be adjusted to optimize titration speed and accuracy.
  - Programmable automatic calculations with individual units
  - Dosing with adjustable dosing and filling rates for optimize speed and accuracy
  - Sample/solution preparation
  - Automatic weight recording when balance is connected
- Scope of delivery:
- Basic unit with magnetic stirrer, with 20 ml exchangeable head, with keyboard, hand controller, stand rod and titration clamp Z 305, power supply 100-240 V