

## TECHNICAL DATA SHEET

---

<b>Part Number</b>	780531
<b>Description</b>	Ebi25te Wireless Temperature Logger With External Sensor -40 To 85c



### Benefits :

- ✓ Wireless temperature recorder, with insert probe ü Measurement range: -40°C +85°C
- ✓ Accuracy: temperature:  $\pm 0.5^{\circ}\text{C}$  (-20°C +40°C)  $\pm 0.8^{\circ}\text{C}$  (-30°C -20°C; +40 °C +60°C)  $\pm 1.5^{\circ}\text{C}$  (-40°C -30°C; +60°C +85°C)
- ✓ Resolution: 0.1°C ü Memory capacity: 288 measured values (per channel)
- ✓ Probe: NTC
- ✓ Measurement interval: 1 min. at 24 hours, adjustable
- ✓ Probe: Ø 3.8 mm ü Cable: Ø 5 mm, L = 50 mm, with 3m PTFE cable
- ✓ Radio frequency: 868 MHz
- ✓ Protection class: IP67
- ✓ Battery: 3.6 V lithium (replaceable by the user)
- ✓ Factory calibration certificate
- ✓ Weight: approx. 65g
- ✓ Dimensions: 95 x 48 x 27 mm

### Sales unit :

1 Ebi25te Wireless Temperature Logger

---