


PDS No. 78109x	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	384 Well Microplate, PS, $\mu$ Clear <sup>®</sup>				 <b>greiner bio-one</b>
	Greiner Item-No. 78109x				
Valid for Item-No.:	<b>781090</b> (sterile)	<b>781091</b> (sterile)	<b>781092</b> (sterile)	<b>781093</b> (sterile)	<b>781094</b> (sterile)
	<b>781095</b>	<b>781096</b>	<b>781097</b> (sterile)	<b>781098</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Microplate, 384 well, $\mu$ Clear <sup>®</sup> , F-bottom (flat), rounded square well design, alphanumeric well coding 781090, -091, -098: physical surface treatment, with lid (single position, low profile), sterile 781092, -093: physical surface treatment, sterile 781094, -097: high binding, sterile 781095, -096: medium binding
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	Total volume: 131 $\mu$ l Working volume: 15 -110 $\mu$ l 781090, -091, -092, -093, -098: growth area: 10 mm <sup>2</sup>
1.4	Material / Resin	Plate: PS (Polystyrene), free of heavy metal Foil: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate: 781090, -091, -092, -096, -097: plate black, bottom clear 781093, -094, -095, -098: plate white, bottom clear Lid: white
1.6	Sterilization	781090, -091, -092, -093, -094, -097, -098: SAL 10 <sup>-3</sup> 781095, -096: no
1.7	Quality Control	- Raw Material-Control: physical and immunological testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - 781090, -091, -098: the covered microplates can be easily gripped on the side in automated processes

2.	Features	
2.1	Basic features	All Item-No.: free of detectable DNase/RNase, human DNA and pyrogens. 781090, -091, -092, -093, -098: contents non-cytotoxic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	781090, -091, -092, -096, -097: 3000 x g: swinging-bucket rotor 781093, -094, -095, -098: 4000 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	4 years after month of production
2.7	Other Information	-

3.	Packaging	781090	781091, -098	781092, -093, -094, -095, -096, -097
3.1	Pieces / Bag	20	8	10
3.2	Pieces / Box	120	32	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality	Certificate of Quality	Certificate of Quality

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 2 December 2014	Date 3 December 2014	Date 3 December 2014	
Date 28.04.2014	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.