

TECHNICAL DATA SHEET

Part Number	783704
Description	1.0ml Tri-coded Tube, 96-format, External Thread



Benefits :

- ✓ Max. working volume :1.0ml
- ✓ Tube height :46.2mm
- ✓ Tube height with cap : 49.6mm
- ✓ Inner diameter : 6.5mm
- ✓ Outer diameter with cap : 8.7mm
- ✓ Center to center : 9.0mm
- ✓ Minimum temperature : -196°C
- ✓ Tube height in rack : 50.5mm
- ✓ Overall rack height (including lid) : 53.2mm
- ✓ Next Generation manufacturing techniques
- ✓ Code readable without removing tubes from racks
- ✓ Suitable for cryogenic storage
- ✓ Made from high-quality virgin PP
- ✓ Sealable with screw caps/TPE septum caps

Sales unit :

960 Tubes