

PRODUCT SPECIFICATION

Product Code:	FR4LF
Product Description:	Flasks, Reaction, Wide Neck, Flat Flange
Brand:	Pyrex
Type of Glass:	Borosilicate Glass Type 3.3, which complies with ISO3585 European and US Pharmacopoeia Glass Type 1
Product Specification:	Nom. Capacity (ml) 4000 Flange Bore (mm) 100 Body O.D. (mm) 207 Height (mm) 225
Pack Quantity:	1 piece per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed: M.Ashley
Quality Technician

Date: 08 Oct 2019