

SteriPlast® sample scoop

The SteriPlast® sample scoop made of polystyrene is produced and assembled in our clean room of class 7 (10,000, except for the 1000 ml scoop), individually packaged for disposable use and then sterilised by gamma rays.

The sample scoop is ideal for powders, granulates and pastes.

The scoop with transparent cover can be closed immediately after sampling and the sterility and exclusion of contamination can therefore be safeguarded.

- ▶ PS, white
- ▶ Individually packaged and sterilised by gamma rays
- ▶ Production, assembly and packaging according to clean room class 7
- ▶ With EU foodstuffs and FDA approval



Ergonomics were particularly important to us for the design of the sample scoop. The scoop stands horizontally due to the special handle design and the slightly angled front edge and can therefore be put down without the sample material trickling out.

SteriPlast® sample scoop									
Capacity ml	Model	Length mm	Pieces per pack	Price €	Quantity scale from	Scale price €	Item no.		
25	without cover	141	10		10		5378-1001	N4	
25	with cover	142	10		10		5378-1002	N4	
50	without cover	160	10		10		5378-1003	N4	
50	with cover	161	10		10		5378-1004	N4	
100	without cover	205	10		10		5378-1005	N4	
100	with cover	206	10		10		5378-1006	N4	
150	without cover	216	10		10		5378-1007	N4	
150	with cover	218	10		10		5378-1008	N4	
250	without cover	232	10		10		5378-1013	N4	
250	with cover	235	10		10		5378-1014	N4	
500	without cover	280	10		10		5378-1015	N4	
500	with cover	283	10		10		5378-1016	N4	
1000	without cover	332	10		10		5378-1019	N4	