

5334-0xxx Volumetric spoons

Food contact compliance statement:

Material complies with EU regulation No. 10/2011 as amended, Annex 1 (Union list of authorized substances)

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up. The supplier – Bürkle GmbH in Bad Bellingen/Germany – is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of certificate is 2284-08.

16. July 2013



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



Material certificate

The media-contacting parts of the articles on the following pages are manufactured from stainless steel or Aluminum, some in combination with PTFE. They are all designed to have a smooth finish and to be internally crevice free. They can be easily cleaned in line with GMP.

The declaration for stainless steel is based on REGULATION (EC) No 1935/2004 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. "The principle underlying this regulation is that any material or article intended to come into contact, directly or indirectly, with food must be sufficiently inert to preclude substances from being transferred to food in quantities large enough to endanger human health or to bring about an unacceptable change in the composition of the food or a deterioration in its organoleptic properties."

Stainless steel is in general widely used and suitable for most foodstuffs. If using highly corrosive or aggressive materials, please check our List of "Chemical-Resistance":

www.buerkle.de/en/chemical-resistance or contact us for suitability.

The declaration is based on our current state of knowledge and information provided by our supplier at the time that the document was drawn up.

The supplier - Bürkle GmbH in Bad Bellingen/Germany - is certified according to the standard DIN EN ISO 9001 by the DQS (German Society for Quality Assurance) since 1995. The number of the certificate is 002284-QM15.

05. July 2023



Bürkle GmbH, Bad Bellingen,
Martin Saint-Denis, Managing Director



[illegible]

[illegible]

[illegible]

Article	Description	Alu	PTFE	PA	1.4034	1.4301 V2A	1.4305	1.4310	1.4571 V4A	1.4404 V4A
5323-0190	Ice borer, V4A, inner-xouter-Ø 8x10mm, depth 10cm								x	x
5323-0200	Ice borer, V4A, inner-xouter-Ø 16x21mm, depth 20cm								x	x
5323-2010	Ice borer, V4A, inner-xouter-Ø 50x54mm, depth 20cm								x	x
5324-0001	PharmaSpoon, V4A, length 650 mm, hand polished									x
5324-0010	PharmaScoop, V4A, 10 ml vol., Ø 25 mm, length 45mm									x
5324-0011	PharmaScoop, V4A, 50 ml vol., Ø 38 mm, length 70mm									x
5324-0012	PharmaScoop, V4A, 100ml vol., Ø 50mm, length 100mm									x
5324-0013	PharmaScoop, V4A, 750ml vol., Ø 95mm, length 210mm									x
5324-0014	PharmaScoop, V4A, 1000 ml vol., Ø110 mm, 235 mm									x
5324-0015	PharmaScoop, V4A, 1250 ml vol., Ø125 mm, 195 mm									x
5324-0016	PharmaScoop, V4A, 1700 ml vol., Ø135 mm, 230 mm									x
5324-0017	PharmaScoop, V4A, 2500 ml vol., Ø185 mm, 215 mm									x
5324-0018	PharmaScoop, V4A, 500ml vol., Ø 95mm, length 140mm									x
5324-0019	PharmaScoop with long handle, V4A, 1000 ml									x
5326-0001	Liquid-Sampler, open with thumb, V4A, 50 ml									x
5326-0002	Liquid-Sampler, open with thumb, V4A, 100 ml									x
5326-0003	Liquid-Sampler, place onto bottle, V4A, 50 ml									x
5326-0004	Liquid-Sampler, place onto bottle, V4A, 100 ml									x
5334-8000	Stainless steel measuring spoons set (1,25-15 ml)					x				
5337-2020	Sampling drill for stalk material, Ø25mm, 200mm				x	x				x
5337-2040	Sampling drill for stalk material, Ø25mm, 400mm				x	x				x
5337-4020	Sampling drill for stalk material, Ø45mm, 200mm	x			x	x				
5337-4040	Sampling drill for stalk material, Ø45mm, 400mm	x			x	x				
5338-0410	StickProof, V4A, lentgth 410 mm, Ø 25mm, ISTA-Norm									x
5338-1100	StickProof, V4A, lentgth 1100mm, Ø 25mm, ISTA-Norm									x
5339-0300	IceSampler, V4A, lentgth 300mm, vol. 50ml, Ø 38mm									x
5339-0550	IceSampler, V4A, lentgth 550mm, vol. 50ml, Ø 38mm									x
5339-1050	IceSampler, V4A, length 1050mm, vol. 50ml, Ø 38mm									x
5340-1100	SiloPicker, V4A, vol. 100 ml, 130 cm, outer-Ø 50mm									x
5341-1002	Suction head, suction drill 2mm, Mini ViscoSampler									x
5341-1004	Suction head, suction drill 4mm, Mini ViscoSampler									x
5341-1006	Suction head, suction drill 6mm, Mini ViscoSampler									x
5341-1060	Mini ViscoSampler, V4A, 60 cm, vol. 30 ml, Ø 15 mm		x							x
5341-1100	Mini ViscoSampler, V4A, 100 cm, vol. 50 ml, Ø 15mm		x							x
5341-1150	Mini ViscoSampler, V4A, 150 cm, vol. 75 ml, Ø 15mm		x							x
5346-0010	PharmaPicker case, complete in transport box									x
5346-1000	PharmaPicker, V4A, w/collection cylinder 0,1-1,2ml									x
5346-1001	Collection cylinder for 0,1-1,2ml, w/o head									x
5346-1010	Volume head, capacity 0,1 ml, PharmaPicker									x
5346-1020	Volume head, capacity 0,2 ml, PharmaPicker									x
5346-1060	Volume head, capacity 0,6 ml, PharmaPicker									x
5346-1100	Volume head, capacity 1 ml, PharmaPicker									x
5346-2000	PharmaPicker, V4A, w/coll. cylinder for 1,25-2,5ml									x
5346-2001	Collection cylinder for 1,25-2,5ml, w/o head									x
5346-2200	Volume head, capcity 2 ml, PharmaPicker									x
5346-2250	Volume head, capcity 2,5 ml, PharmaPicker									x
5346-3010	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,01ml								x	x
5346-3050	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,05ml								x	x

Article	Description	Alu	PTFE	PA	1.4034	1.4301 V2A	1.4305	1.4310	1.4571 V4A	1.4404 V4A
5346-3100	Powder sampler Mini, AISI316L, Ø6mm, 300mm, 0,1ml								x	x
5347-0005	Scratcher, V4A, 5 ml, flute length 22 mm, Ø 25 mm									x
5347-0010	Scratcher, V4A, 10 ml, flute length 45 mm, Ø 25 mm									x
5347-0025	Scratcher, V4A, 25 ml, flute length 110 mm, Ø 25mm									x
5351-1000	All-layer-scratcher, V4A, Ø 25mm, 100cm, vol.180ml		x			x				x
5354-1000	Stainless steel scoop, V2A, 115cm, capacity 1000ml					x				
5354-2000	Stainless steel beaker, V2A, 1000 ml, TeleScoop					x				
5354-2010	Stainless steel beaker, extended, 1000ml, TeleScoop					x				
5354-3000	Pendulum beaker V2A, capacity 1000 ml, TeleScoop					x	x			
5356-1500	EasyScoop, aluminium, 1500 ml, inner grip	x								
5356-2000	EasyScoop, aluminium, 2000 ml, inner grip	x								
5357-0001	Sampler Liquid CupSampler, V4A, 100 ml,100 cm									x
5358-0000	Sample insert, 0,0 ml, Novartos Multi									x
5358-0005	Sample insert, 0,5 ml, Novartos Multi									x
5358-0010	Sample insert, 1 ml, Novartos Multi									x
5358-0015	Sample insert, 1,5 ml, Novartos Multi									x
5358-0020	Sample insert, 2 ml, Novartos Multi									x
5358-0100	Sample insert, 0,0 ml, Novartos Uno									x
5358-0105	Sample insert, 0,5 ml, Novartos Uno									x
5358-0110	Sample insert, 1 ml, Novartos Uno									x
5358-0120	Sample insert, 2 ml, Novartos Uno									x
5358-0130	Sample insert, 3 ml, Novartos Uno									x
5358-0140	Sample insert, 4 ml, Novartos Uno									x
5358-0600	Novartos Uno, V4A, 600 mm, 2 chambers, w/o inserts									x
5358-0602	Novartos Multi, 60cm,double,8 chambers,w/o inserts									x
5358-0603	Novartos Multi,60cm,triple,12 chambers,w/o inserts									x
5358-0630	Novartos Uno, V4A, 600 mm, 3 chambers, w/o inserts									x
5358-1000	Novartos Multi, V4A, 100cm, 14 chambers, w/o inserts									x
5358-1001	Novartos Uno, V4A, 1000 mm, 2 chambers, w/o inserts									x
5358-1003	Novartos Multi,100cm,triple,21 chambers,w/o inser.									x
5358-1011	Novartos V4A, head									x
5358-1014	Spitze, V4A, abgerundet R 5 mm, Novartos									x
5358-1030	Novartos Uno, V4A, 1000 mm, 3 chambers, w/o inserts									x
5358-1400	Novartos Uno, V4A, 1400 mm, 2 chambers, w/o inserts									x
5358-1402	Novartos Multi,140cm,double,18 chambers,w/o inser.									x
5358-1403	Novartos Multi,140cm,triple,27 chambers,w/o inser.									x
5358-1430	Novartos Uno, V4A, 1400 mm, 3 chambers, w/o inserts									x
5359-1600	Valve lance, V4A, length 550 mm, Ø 32 mm									x
5359-1800	Valve lance, V4A, length 800 mm, Ø 32 mm									x
5360-0600	EasySampler V4A 60ml									x
5360-1000	EasySampler, V4A, length 104,5 cm, volume 200 ml									x
5369-0300	Pigment lance, flute open.25%, length 300mm, 145ml									x
5369-0310	Pigment lance, flute open. 50%, flute 300mm, 85 ml									x
5369-0600	Pigment lance, flute open. 25%, flute 600mm, 312ml									x
5369-0610	Pigment lance, flute open. 50%, flute 600mm, 179ml									x
5369-0900	Pigment lance, flute open. 25%, flute 900mm, 480ml									x
5369-0910	Pigment lance, flute open. 50%, flute 900mm, 273ml									x
5370-0120	Scoop AISI 304, V2A, capacity 120 ml, Ø 50 mm					x				

Article	Description	Alu	PTFE	PA	1.4034	1.4301 V2A	1.4305	1.4310	1.4571 V4A	1.4404 V4A
5370-0450	Scoop AISI 304, V2A, capacity 450 ml, Ø 75 mm					x				
5370-1000	Scoop AISI 304, V2A, capacity 1000 ml, Ø 100 mm					x				
5370-1500	Scoop AISI 304, V2A, capacity 1500 ml, Ø 120 mm					x				
5370-2000	Scoop AISI 304, V2A, capacity 2500 ml, Ø 150 mm					x				
5371-0100	FoodScoop, V2A, 100ml, total l:215mm, LxW 120x60mm					x				
5371-0200	FoodScoop, V2A, 200ml, total l:230mm, LxW 140x75mm					x				
5371-0250	FoodScoop, V2A, 250ml, total l:245mm, LxW 155x85mm					x				
5371-0350	FoodScoop, V2A, 350ml, total: 290mm, LxW 180x100mm					x				
5371-0500	FoodScoop, V2A, 500ml, total: 350mm, LxW 220x130mm					x				
5380-0506	Extension rod, 500 mm, M6 thread, PharmaPicker									x
5380-0508	SiloPicker extension 50 cm									x
5380-1006	Extension rod, 1000 mm, M6 thread, PharmaPicker									x
5380-1008	SiloPicker extension 100 cm									x
5382-0040	Sampler Tubus, V4A, 40 cm, Ø 40 mm, ISTA standard									x
5382-0080	Sampler Tubus, V4A, 80 cm, Ø 40 mm, ISTA standard									x
5383-1000	Spatula stainless steel, V4A, 300 mm									x
5386-0001	Micro spatula stainless steel, V2A, length 170 mm					x				
5386-0010	Palette knife spatula stainl. steel, LxW 200x21 mm					x				
5386-0012	Palette knife spatula stainl. steel, LxW 250x25 mm					x				
5386-0022	Spoon spatula stainless steel, length 180 mm, 4 ml					x				
5386-0025	Spoon spatula stainless steel, length 300 mm, 11ml					x				
5386-0030	Scraper stainless steel, AISI 316, LxW 230 x 75mm									x
5386-0035	Scraper stainless steel, AISI 316, LxW 270 x 100mm									x
5386-0100	Forceps, stainless steel, sharp, bent form, 105 mm					x				
5386-0102	Forceps, stainless steel, sharp, bent form, 130 mm					x				
5386-0104	Forceps, stainless steel, sharp, bent form, 160 mm					x				
5386-0200	Forceps, stainless steel, sharp, straight, 105 mm					x				
5386-0202	Forceps, stainless steel, sharp, straight, 130 mm					x				
5386-0204	Forceps, stainless steel, sharp, straight, 160 mm					x				
5386-0300	Forceps, stainless steel, blunt, straight, 105 mm					x				
5386-0302	Forceps, stainless steel, blunt, straight, 130 mm					x				
5386-0304	Forceps, stainless steel, blunt, straight, 160 mm					x				
5388-1000	BeefSteaker, V2A, Ø 25mm				x					
5388-1010	Drill head for carcass sampling				x					
9604-0001	Funnel, stainless steel, V2A, outer-Ø 60 mm					x				
9604-0003	Funnel, stainless steel, V2A, outer-Ø 100 mm					x				
9604-0005	Funnel, stainless steel, V2A, outer-Ø 150 mm					x				
9604-0008	Funnel, stainless steel, V2A, outer-Ø 300 mm					x				