

TECHNICAL DATA SHEET

Part Number	79606-500ML-HON
Description	Phosphoric acid

Technical specifications :

- ✓ phosphoric acid, for hplc, 85-90% 500ml
- ✓ CAS N°7664-38-2
- ✓ Chemical purity 85-90%
- ✓ Molecular weight 98 g/mol
- ✓ Molecular formula H₃O₄P
- ✓ Density 1.71 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 500mL

- ✓ GHS pictograms GHS05
- ✓ H Statements H290 ; H314
- ✓ P Statements P284 ; P280 ; P260 ; P234

Sales unit :

500mL Phosphoric acid
