

Technical Datasheet

Article no.: 797

Tetrahydrofuran for HPLC (min. 99.9 % (unstab.))

C₄H₈O

For laboratory use.

CAS-No: 109-99-9
Appearance/condition: Clear, colourless liquid
Melting point: -108 °C
Boiling point: 64 – 66 °C
Density (20 °C): 0.89 g/ml
Molar mass: 72.11 g/mol

Parameter	Value
Refractive index (20 °C)	1.405 – 1.409
Water (KF)	max. 200 mg/kg
Non-volatile substances	max. 5 mg/kg
Free acids (as CH ₃ COOH)	max. 20 mg/kg
Peroxides (as H ₂ O ₂)	max. 300 mg/kg
UV transmittance at 240 nm	min. 20.0 %
UV transmittance at 250 nm	min. 45.0 %
UV transmittance at 300 nm	min. 90.0 %
UV transmittance at 320 nm	min. 95.0 %
Colour (APHA)	max. 10
Filtered through 0.2 µm.	

Version No. 02/2024-10-14