

Product Specification

Product Code	CP0050P
Product Description	AZLON MEASURING CYLINDER, CLASS B, MOULDED GRADUATIONS, PP, 50 ML
Material particular	Polypropylene
Capacity ml	50
Colour	Ultra Clear PP
Autoclavable	Yes, See www.scilabware.com
Pack quantity	10

Dimensions and tolerances

Height mm	200
Diameter mm	26
Graduation ml	10.0

Standards

Conforms to	ISO 6706
-------------	----------

Product Compliance

TSE/BSE Compliant	Yes
Reach	Yes
Latex	Yes

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO 9001:2008

Edel Tuite

Azlon Quality Manager

Date issued 1/19/2017