

Product Specification

| | |
|---------------------|---|
| Product Code | CP0250P |
| Product Description | AZLON MEASURING CYLINDER, CLASS B, MOULDED GRADUATIONS, PP, 250 ML |
| Material particular | Polypropylene |
| Capacity ml | 250 |
| Colour | Ultra Clear PP |
| Autoclavable | Yes, See www.scilabware.com |
| Pack quantity | 1 |

Dimensions and tolerances

| | |
|---------------|------|
| Height mm | 315 |
| Diameter mm | 42 |
| Graduation ml | 20.0 |

Standards

| | |
|-------------|----------|
| Conforms to | ISO 6706 |
|-------------|----------|

Product Compliance

| | |
|-------------------|-----|
| TSE/BSE Compliant | Yes |
| Reach | Yes |
| Latex | Yes |

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO 9001:2008

Edel Tuite

Azlon Quality Manager

Date issued 1/19/2017