

Soft containers



- Material: Nylon 210/210, PVC cristal transparent
- Size 30 x 9 x 16 cm
- Various colors available
- Zipper closure
- Three faces in transparent PVC
- A face with Velcro

REF	Color
BS2405A	Orange
BS2405B	Blue
BS2405G	Yellow
BS2405R	Red
BS2405V	Green



Biosigma s.r.l.

a dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |
Fax ++39 0426 748030 | Email info@biosigmaeu.com | www.biosigma.it