

ATTESTATION OF CONFORMITY

Issued to: Hangzhou Bioer Technology Co., Ltd.
1192 Bin'an Road, Hi-tech (Binjiang) District, Hangzhou, Zhejiang, 310053 P.R.
China

For the product: Fluorescent Quantitative PCR Detection system

Trade name: BIOER

Type/Model: FQD-96C

Ratings: 200— 240V~, 50/60 Hz, 1000 W

Manufactured by: Hangzhou Bioer Technology Co., Ltd.
1192 Bin'an Road, Hi-tech (Binjiang) District, Hangzhou, Zhejiang, 310053 P.R.
China

Requirements: EN 61326-1:2013;
EN IEC 61000-3-2:2019
EN 61000-3-3:2013+A1:2019

This Attestation is granted on account of an examination by DEKRA, the results of which are laid down in a confidential file no. 6102431.50

This Attestation implies that the examined types are in accordance with the standards designated under the Electromagnetic compatibility directive 2014/30/EU.

The examination has been carried out on one single specimen or several specimens of the product, submitted by the manufacturer. The Attestation does not include an assessment of the manufacturer's production. Conformity of his production with the specimen tested by DEKRA is not the responsibility of DEKRA.

The CE marking may be affixed on the product if all relevant and effective EC directives are complied with.

Arnhem, 18 May 2021

Number: 6102431.01AOC

DEKRA Testing and Certification (Shanghai) Ltd.



Kreny Lin
Certification Manager

© Integral publication of this attestation and adjoining reports is allowed

