



Item Number / Product Name / Quantity
Country of Origin: Taiwan (ROC)

6803	5-Hole Multi-Purpose Rack	100 units / carton
-------------	----------------------------------	---------------------------

Inspection Item	Specification	Testing Method	Result	Judgment
Appearance	Must be consistent 5	Visual Inspection 6	Consistent 7	OK
Color	Match color card 8	Visual Inspection 9	Match 10	OK
Length	101 mm 11	Vernier Caliper 12	+/- 0.3 mm 13	OK
Width	50 mm 14	Vernier Caliper 15	+/- 0.3 mm 16	OK
Height	96 mm 17	Vernier Caliper 18	+/- 0.3 mm 19	OK
Hole Diameter	29.5 mm 20	Vernier Caliper 21	+/- 0.3 mm 22	OK
Cracks/Damage	No cracks or damage 23	Visual Inspection 24	None 25	OK
Foreign Matter	No obvious foreign matter 26	Visual Inspection 27	None 28	OK