




Slip and lock the racks together.



Stackable.


Two clear lids included.



5405-59



Revolve to create the perfect workstation.



5305-59

Locks in place

### Sliplock Flipper

- Expandable mini workstation, includes 4 interchangeable racks and 2 clear lids.
- Capacity:
  - 16 x 0.5 mL tubes per rack.
  - 24 x 1.5/2.0 mL tubes per rack.
  - 60 x 0.1/0.2 mL PCR tubes per rack.
  - 12 x 8-strip or 11 x 12-strip PCR tubes per 2 racks.
  - 1 x unskirted or semi-skirted 96-well PCR plate per 2 racks.
- Dimensions (L x W x H):
  - 1 rack: 134 x 54 x 35mm
  - 2 racks: 134 x 108 x 35mm (134 x 110 x 58 with lid).
  - 4 racks: 215 x 134 x 35mm (220 x 134 x 58 with lids).
- Temperature range: -80°C to +121°C.
- Rugged polypropylene construction.

Item #	Description	Colors	Packaging
5305-59	Sliplock Flipper 4 Racks and 2 Clear Lids	Green, Blue, Red, Yellow	1 workstation/unit, 4 units/case

### Revolver Flipper - the ultimate workstation!

- Add or subtract rack units to suit your needs.
- PCR workstation options:
  - 48 x 0.2 mL tubes per block.
  - 6 x 8 tube strips per block.
  - 8 x 12 tube strips per 2 blocks.
  - 1 x non-skirted or semi-skirted 96 well plate per 2 blocks.
- Four rotate-and-lock racks.
- Fast reconfiguration for different tube sizes and types.
- Temp range: -80°C to +121°C.
- Rugged polypropylene construction.

Item #	Dimensions	Tube Options	Colors	Packaging
5405-59	Revolver Flipper 257 x 108 x 108 mm	Side 1: 6 x 0.5 mL and 12 x 1.5/2.0 mL Side 2: 8 x 5 or 15 mL Side 3: 2 x 50 mL Side 4: 48 x 0.2 mL	Green, Blue, Red, Yellow	4 rack blocks/unit, 4 units/case