

## 3247-00 0.2 mL 8-Strip PCR UltraFlux® i Tubes, Flat Cap



### 0.2 mL 8-Strip PCR UltraFlux® i Tubes, Flat Cap

**SKU:** 3247-00

**Cap Description:** Flat, Optically-Clear, Angularly Attached

**Tube Description:** Strip of 8

**Tube Colors:** -00 Clear, -01 Blue, -02 Green, -04 Red, -05 Violet, -06 Yellow, -09 Assorted, -40 White

**Packaging:** 120 strips/unit, 10 units/case

3247-00, 3247-00S, 3247-01, 3247-01S, 3247-02, 3247-02S, 3247-04, 3247-04S, 3247-05, 3247-05S, 3247-06, 3247-06S, 3247-09, 3247-09S, 3247-40, 3247-40S

- Robust inter-well linkages make the strips rigid and easy to handle, preventing spills.
- Angularly-attached caps prevent lids and hinges from interfering with each other, allowing for easier operation in an SBS format rack.
- Available with clear, white or colored tubes.
  - For RT-PCR, choose the 3247-40 with white tubes and optically-clear flat clear caps. A unique two-step process creates the perfect environment for consistent, reliable qPCR performance. White wells offer enhanced fluorescence reflection and optically-clear caps ensure that the fluorescence signal is uninhibited.
  - Also available in assorted colors for easy organization.
- Manufactured from high-quality, virgin polypropylene to ensure purity and provide excellent optical clarity.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.