

## Block as tempering tub with lid

(Article number 800012100)

for vials, flasks, deep well plates or MTPs

Exchangeable thermoblock as tempering tub for Deep-Well-plates or other tubes, 128 x 85 x 70 mm.

### Key benefits:

1. Safety: All blocks are autoclavable
2. Stability: Strong block construction
3. Plug and Play: simple and fast exchanging - blocks suitable for all types of devices
4. Various applications: Accessories like Insulating top, Anti-condensation-plate and much more

### Photos



### Technical Details

#### Technical Data

#### BW 01 (Art. 800012100)

Dimension metal block 128,5 x 85,5 x 25 mm

Dimension plastic insulation 145 x 102 x 98 mm

Tubes	universal, vials, flasks, deep well plates or MTPs
Number of drillings	none
Drilling diameter	128,5 x 70mm
Depth of drillings	73 mm
Shape of drillings	flat
Weight	0,9 kg

## **Accessories**

### **Suitable accessories for this block:**

Tempering cover for systems with one block (Article number 800013800)

Tempering cover for systems with two blocks (Article number 800013900)

Insulating lid is already included