



8 l / min (air)

700 mbar abs.



Best price-performance ratio

This very quiet membrane pump generates a final vacuum of 700 mbar and aspirates 8 l / minute. The Aspiration System AZ is perfectly suited to aspirate fast small amounts of liquid. It is particularly well suited to aspirate fluids from floating cells. The bottleholder and footswitch are an option for comfortable handling.

The security filter prevents the leaking of liquid.

Polypropylen bottles for 2 or 4 litres with quick fasteners and ventilation applications for autoclaving are available.

Components:

- » Membrane vakuumpumpe (8 l / min air, 700 mbar abs.)
- » 2l- or 4l-bottle incl. fast lock clutches, filter, ventilation insert and 3mm silicone hose

Dimensions: 75 x 170 x 60 mm

Weight: 0,8 kg

Electric power output: 5 W

Power supply: 230V, 50Hz / 115V, 50...60 Hz

Article no.	Article description
981523101	AZ 02 Aspiration system with 2 l bottle
981523001	AZ 04 Aspiration system with 4 l bottle

Key benefits:

1. High aspiration volume
2. Extreme small instrument dimensions