



Enjoy the benefits without the pump

When a vacuum port is already installed in the lab, the Aspiration Systems AE and AK are suitable for connection to the existing pump.

The security filter prevents the leaking of liquid.

Polypropylen bottles for 2 or 4 litres with quick fasteners and ventilation applications for autoclaving are available.



Components AE 02 / 04

- » 2l- or 4l-bottle incl. fast lock clutches, filter, ventilation insert and 3mm silicone hose

Components AK 02 / 04

- » 2l- or 4l-bottle incl. fast lock clutches, filter, ventilation insert and 3mm silicone hose
- » bottleholder
- » handle

Article no.	Article description
900015400	AE 02 Aspiration system with 2 l bottle
900015500	AE 04 Aspiration system with 4 l bottle
900014800	AK 02 Aspiration system with 2 l bottle
900014900	AK 04 Aspiration system with 4 l bottle