

Description

- PLASTIC DISINFECTION MAT 90x60x1.2cm

Part number

- 853018

Pieces by sales : 1



Description

- Disinfection mat, effectively stop transported contamination
- Pour our disinfectant, after diluting it, onto the carpet.
- A simple pass over this mat disinfects cart wheels and shoe soles.
- Due to its effectiveness in countering transported contamination, this mat is increasingly used in collective and commercial catering, in the food industry, in the catering trade, in hotels, restaurants, laboratories...
- The mat is placed at the entrance to production areas, workshops, etc. to maintain cleanliness and reduce
- Risk of contamination. Just put it down. It doesn't slip.
- Features
 - sizes : Rug 90 x 60 x 1.2 cm: ref. 002016.
 - 100% plastic mat, which makes it easier to maintain.
- Mat with 4 beveled sides, which makes it easy to get on and off the mat.
- Manual
 - Carpet 90 x 60 x 1.2 cm: dilute 12 ml of disinfectant in 2 to 2.4 L of water then pour the mixture onto the carpet.
 - Refill regularly with disinfectant depending on traffic.
- RECOMMENDATIONS
 - When receiving the rug, please do not open the package with a cutter. Do not fold the mat to dry.

