

ClearLine stomacher bag

Description :

- ✓ Capacity : 400 ml
- ✓ Dimensions W x H : 180 x 300 mm
- ✓ Matériau : polyéthylène
- ✓ Sterile
- ✓ Without filter
- ✓ Flexible
- ✓ High mechanical strength
- ✓ Storage temperature : 15 to 30°C
- ✓ Can be used at 37°C for a maximum of 24 hours
- ✓ Mixing mustn't exceed 3min with a overhead stirrer
- ✓ Not suitable for very pasty products

Unités :

500 bags