

TECHNICAL DATA SHEET

Part Number	85815-1KG-HON
Description	Silica gel Rubin

Technical specifications :

- ✓ silica gel rubin, drying agent, with moisture indicator, cobalt free 1kg
- ✓ CAS N°112926-00-8
- ✓ Density 2 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form granules

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 1kg

- ✓ P Statements P280

Sales unit :

1kg Silica gel Rubin
