

**COMMERCIAL CENTER&R-D Center**  
 #203, Building 2, 89 Shengli Rd, Jiangning  
 District, Nanjing 211106, China  
 Tel.: 0086 25 86216803  
 Fax: 0086 25 86214463  
 E-mail: info@citotest.com

**PRODUCTION**  
 #339, Beihai Road, Haimen City,  
 Jiangsu 226100, P.R.China  
 Tel.: 0086 513 82259348  
 Fax: 0086 513 82110195  
 E-mail: info@citotest.com.cn



**CITOTEST LABWARE MANUFACTURING CO.,LTD** [www.citotest.com](http://www.citotest.com)

# PRODUCT DATA SHEET


## GENERAL INFORMATION

**PRODUCT NAME:** Pipette Tip  
**FURTHER DESCRIPTION:** PIPETTE TIPS WITH FILTER, NATURAL, PP MATERIAL. VOL. 200µL/10-200µL, MOLDED GRADUATION AT 10/50/100µL.UNIVERSAL FIT. RNASE, DNASE AND PYROGEN FREE. STERILE, 96 TIPS/RACK  
**P/N:** 4330-7144  
**CATEGORY NUMBER:** \_\_\_\_\_


## SPECIFICATION

**VOLUME (ml):** 200µL/20-200µL  
**COLOR:** Natural  
**MATERIAL:** PP

## PRODUCT PICTURE:



## TECHNICAL PROCESS:



## INTENDED APPLICATION

Used to transport a measured volume of liquid.

## PACKING

96 Tips/Rack, 10 racks/inner box, 5x10 racks/carton

## BOX:



50\*25mm

Pipette tips with filter, Natural, PP material,  
 Vol. 200µL, Universal fit, RNase, DNase and  
 Pyrogen Free, Sterile, 96 Tips/Rack


**4330-7144**

**96pcs**

INNERBOX: [ Innerbox Label Size(mm) 100x70 ]



CARTON: [ Carton Label Size(mm) 100\*130mm ]



**STERILIZATION**

- STERILE EO  
  STERILE R  
  STERILE A

**STORAGE INSTRUCTION**

Store in a cool, dry place.

**TERM OF VALIDITY**

60 month

**REMARKS**

DOMINIQUE DUTSCHER SAS