

COMMERCIAL CENTER&R-D Center
 #203, Building 2, 89 Shengli Rd, Jiangning
 District, Nanjing 211106, China
 Tel.: 0086 25 86216803
 Fax: 0086 25 86214463
 E-mail: info@citotest.com

PRODUCTION
 #339, Beihai Road, Haimen City,
 Jiangsu 226100, P.R.China
 Tel.: 0086 513 82259348
 Fax: 0086 513 82110195
 E-mail: info@citotest.com.cn



CITOTEST LABWARE MANUFACTURING CO.,LTD www.citotest.com

PRODUCT DATA SHEET


GENERAL INFORMATION

PRODUCT NAME: Pipette Tip
FURTHER DESCRIPTION: PIPETTE TIPS WITH FILTER, PP MATERIAL. VOL. 1000µL/100-1000µL, MOLDED GRADUATION AT 50/100/200/500µL.UNIVERSAL FIT. RNASE, DNASE AND PYROGEN FREE. STERILE, 100 Tips/box
P/N: 4330-7152
CATEGORY NUMBER: _____


SPECIFICATION

VOLUME (ml): 1000µL/100-1000µL
COLOR: Natural
MATERIAL: PP

PRODUCT PICTURE:



TECHNICAL PROCESS:



INTENDED APPLICATION

Used to transport a measured volume of liquid.

PACKING

96 Tips/Rack, 10 racks/inner box, 5x10 racks/carton



BOX:

50*25mm


Pipette tips with filter, Natural, PP material,
 Vol. 1000µL, Universal fit, RNase, DNase and
 Pyrogen Free, Sterile, 100 Tips/Box

4330-7152 **100pcs**

INNERBOX: [Innerbox Label Size(mm) 100x70]



CARTON: [Carton Label Size(mm) 100*130mm]



STERILIZATION

STERILE EO
 STERILE R
 STERILE A

STORAGE INSTRUCTION

Store in a cool, dry place.

TERM OF VALIDITY

60 month

REMARKS

DOMINIQUE DUTSCHER SAS