

CITOTEST[®]

Technical Data Sheet

First edition

Document Number: TDS100127

Cryotubes

1. Technical information	
Product name	Cryotube (PP material)
2. Product range	
CAT. No.	Description
4004-1725/4004-1705/4004-1744/4004-1706/4004-1745/4004-1736	<ul style="list-style-type: none"> ● Vol. 2.0ml ● DIM. (OD X H): 48.5 x12.4mm ● Graduations at 0.5/1.0/1.8ml
4004-1726/4004-1707/4004-1746/4004-1708/4004-1737	<ul style="list-style-type: none"> ● Vol. 5.0ml ● DIM. (OD X H): 91.3 x12.4mm ● Graduations at 1.0/2.0/3.0/4.0/4.5ml
4004-1721/4004-1714/4004-1715	<ul style="list-style-type: none"> ● Vol. 1.0ml ● DIM. (OD X H): 41 x12.3mm ● Graduations at 0.5/1.0ml
4004-1722/4004-1741/4004-1716/4004-1742/4004-1717	<ul style="list-style-type: none"> ● Vol. 2.0ml ● DIM. (OD X H): 49 x12.3mm ● Graduations at 0.5/1.0/1.8ml
4004-1723/4004-1718/4004-1719	<ul style="list-style-type: none"> ● Vol. 4.0ml ● DIM. (OD X H): 70.3 x12.3mm ● Graduations at 1.0/2.0/3.0/3.6ml
4004-1738/4004-1739/4004-1740	<ul style="list-style-type: none"> ● Vol. 5.0ml ● DIM. (OD X H): 91.3 x12.4mm ● Graduations at 1.0/2.0/3.0/4.0/4.5ml
3. Intended Use	
CITOTEST Cryotube made of PP material are designed for the storage or/and transportation of biological specimens.	
4. Features	

- Leak-proof, silicon liner mounted in cap
- Internal thread tubes can minimize liquid retention and suitable for vapor phase storing.
- Minimum temperature: -196°C.
- RNase free, DNase free and Non-Pyrogen
- With large white writing surface for specimen identification.
- Autoclavable at 121°C for 15min.

5. Certification

CITOTEST certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2008 Quality Management System.

6. Specification



Type A



Type B



Type C

6.1 External thread



P/N	Type	Packing
4004-1725	A	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1705	A	100pcs/bag, 1000pcs/carton
4004-1744	B	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1706	B	100pcs/bag, 1000pcs/carton
4004-1745	C	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1736	C	100pcs/bag, 1000pcs/carton
4004-1726	A	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton
4004-1707	A	100pcs/bag, 1000pcs/carton
4004-1746	B	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton
4004-1708	B	100pcs/bag, 1000pcs/carton
4004-1737	C	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton

6.2 Internal thread



P/N	Type	Packing
4004-1721	A	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1714	B	100pcs/bag, 1000pcs/carton
4004-1715	C	100pcs/bag, 1000pcs/carton
4004-1722	A	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1741	B	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1716	B	100pcs/bag, 1000pcs/carton
4004-1742	C	50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton
4004-1717	C	100pcs/bag, 1000pcs/carton
4004-1723	A	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton
4004-1718	B	100pcs/bag, 1000pcs/carton
4004-1719	C	100pcs/bag, 1000pcs/carton
4004-1738	A	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton
4004-1739	B	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton
4004-1740	C	25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton

7. Sterilization

STERILE EO
 STERILE R
 STERILE A
 Non-sterile

8. Precautions

- To be used by adequately trained and qualified personnel.
- Not to be used in the liquid phase of liquid nitrogen.
- Disposable of used tubes in suitable disposal unit or in accordance with local regulations.

9. Storage and Shelf Life

Store in a clean, dry environment at room temperature 15-30°C.

10. Version Information

Unless otherwise specified, no part of this document may be reproduced or utilized in any form or by any means, electronic or mechanical, including photocopying and microfilm, without permission from CITOTEST LABWARE MANUFACTURING Co., Ltd.

Citotest Labware Manufacturing Co., Ltd
 339 Beihai Road, Haimen City226100, Jiangsu, P.R. China
 TEL:0086-513-82259348 FAX:0086-513-82110195

Web page: <http://www.citotest.com>