


Test Verification of Conformity

Verification Number: 181200687SHA-V1

On the basis of the referenced test report(s), sample(s) tested of the below product have been found to comply with the standards harmonized with the directives listed on this verification at the time the tests were carried out. Other standards and Directives may be relevant to the product. This verification is part of the full test report(s) and should be read in conjunction with it <them>.

Once compliance with all product relevant  mark directives are verified, including any relevant e.g. risk assessment and production control, the manufacturer may indicate compliance by signing a Declaration of Conformity themselves and applying the mark to products identical to the tested sample(s).

Applicant Name & Address:	Hangzhou Bioer Technology Co., Ltd 1192 Bin'an Road, Hi-tech (Binjiang) District, Hangzhou, Zhejiang, 310053 P.R. China
Manufacturer Site Name & Address	Hangzhou Bioer Technology Co., Ltd 1192 Bin'an Road, Hi-tech (Binjiang) District, Hangzhou, Zhejiang, 310053 P.R. China
Product Description: Ratings & Principle Characteristics:	Fluorescent Quantitative Detection system 100-240V,50/60Hz
Models/Type References:	FQD-96C
Brand Name(s):	BIOER
Standard(s)/Directive(s):	EN61326-1:2013 the EMC directive (2014/30/EU)
Verification Issuing Office Name & Address:	Intertek Testing Services Shanghai Building No.86, 1198 Qinzhou Road (North), Caohejing Development Zone, Shanghai 200233, China
Test Report Number(s):	181200687SHA-001



Signature

Name: Jonny Jing

Position: Manager

Date: Dec 21, 2018