

# 25 mm Microanalysis Glass Filtration Apparatus



VF5



VF8

## Applications


Ideal for small volume vacuum filtration for biological or particulate contamination analysis.

## Specifications

Material	VF5	VF8
Funnel	Borosilicate glass	Borosilicate glass
Support base	Borosilicate glass	Borosilicate glass
Membrane support	Sintered glass	SS screen with PTFE gasket
Clamp	Aluminum	Aluminum
Stopper	Silicone rubber (No.5)	Silicone rubber (No.5)
Receiver flask	Borosilicate glass	Borosilicate glass

## Specifications

Funnel capacity	15 ml	15 ml
Flask capacity	125 ml	125 ml
Filter diameter	25 mm	25 mm
Effective filtration area	2.2 cm <sup>2</sup>	2.2 cm <sup>2</sup>
Hose barb	5/16 inch (8 mm)	5/16 inch (8 mm)



VF5

VF8

## Ordering Information

VF5, 25mm Microanalysis Glass Filtration Set 167200-05

15 ml Glass filter holder (1+2+3) 167120-32

1	Glass funnel, 15 ml	167220-15
2	Aluminum clamp	167240-05
3	Sintered glass support base (Including No.5 silicone stopper)	167230-04

4 Receiver flask, 125 ml 167250-05

VF8, 25 mm Microanalysis Glass Filtration Set 167200-08

15 ml Glass filter holder (1+2+3+4+5) 167120-34

1	Glass funnel, 15 ml	167220-15
2	Aluminum clamp	167240-05
3	SS support screen	167230-32
4	PTFE gasket	167230-42
5	Glass support base (Including No.5 silicone stopper)	167230-08

6 Receiver flask, 125 ml 167250-05