


# UltraSlim UV transilluminator

*As Small As it can be!*

Cat No. UUV-01



-  **Focus on mini-gel:** 150\*80mm viewing
-  **High Performance:** Sensitivity > 1ng
-  **Small and Economic:** Smaller and cheaper




 **MAESTROGEN**

[www.maestrogen.com](http://www.maestrogen.com)  
[info@maestrogen.com](mailto:info@maestrogen.com)



*Detachable UV shield*



**ONE-YEAR  
WARRANTY**

*UV shield for  
fluorescence observation*

**LIGHT, EASY  
TO USE**



<b>Dimension</b>	<b>247 x 132 x 48.5 mm</b>
<b>Wavelength</b>	<b>302 nm</b>
<b>UV Tube</b>	<b>T5 6W UVB</b>
<b>Weight</b>	<b>1.4 Kg</b>
<b>Gel Size</b>	<b>125X60 mm</b>
<b>Power</b>	<b>12V, 1.5A</b>

 **MAESTROGEN**

[www.maestrogen.com](http://www.maestrogen.com)  
[info@maestrogen.com](mailto:info@maestrogen.com)

