

TECHNICAL DATA SHEET

Part Number	89152-250G-HON
Description	Sodium perchlorate monohydrate

Technical specifications :

- ✓ sodium perchlorate monohydrate, for hplc =99,0% 250g
- ✓ CAS N°7791-07-3
- ✓ Chemical purity ≥99%
- ✓ Molecular weight 140.46 g/mol
- ✓ Molecular formula NaClO₄.1H₂O
- ✓ Density 2.02 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form solid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 250g

- ✓ GHS pictograms GHS07 ; GHS08 ; GHS03
- ✓ H Statements H271 ; H302 ; H319 ; H373
- ✓ P Statements P280 ; P260 ; P210

Sales unit :

250g Sodium perchlorate monohydrate
