

## UNITED POLY ENGINEERING PVT. LTD.

D-13/3, OKHLA INDUSTRIAL AREA  
PHASE – II, NEW DELHI, INDIA

[www.unitedpoly.com](http://www.unitedpoly.com)



### SURGICAL BOX, STAINLESS STEEL



#### STAINLESS STEEL (AISI 304 GRADE)

<u>CODE</u>	<u>DIM (MM)</u>
UPL- 214	200 x 90 x 30
UPL -216	200 x 90 x 50

#### TECHNICAL SPECIFICATION:

Basic equipment for nursing and surgical care.

Commonly used for collection, storage, transportation, sterilization of medical/surgical material/instruments/Specimen/Waste etc.

This item can be sterilized in a steam sterilizer or hot air sterilizer if used for sterile purposes

Material – Austenitic Stainless Steel, smooth surface.

Austenitic Stainless Steel composition: approx 8 to 10 % nickel, 18 to 20% chromium.

Thickness: Base 0.8mm, Lid 0.7mm,

Length Tolerance:  $\pm 3$ mm

Cleaning Instructions:

Wash with mild detergent and warm water or place in an autoclave for sterilization.

Do not use abrasive cleaners to maintain the integrity of the surface.