

Product Description

Catalog Number: 3194-335-001-9

Product Description: 50mL Metal-Free Centrifuge Tube
Polypropylene Sterile. Printed Graduations
Flat Natural-Colored Cap

Product Dimensions:

Length of Tube: 4.51" (114.55mm)
Outer Diameter of Tube 1.16" (29.46mm)
Inner Diameter of Tube 1.07" (27.17mm)

Packaging: Racked (cardboard fiber racks)
25 tips per rack (with attached caps)
2 racks per pack (50 tubes)
10 packs per case (500 tubes total)

Autoclavable. Freezable. Tubes can be centrifuged to 20,000 RCF.

MicroBiological testing:

Certified testing is performed lot by lot. These tests include Nuclease (RNase and DNase), Endotoxin (pyrogen) and ATP. Certificates are on file at the production facility.

Performance Testing:

Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.

Sterility:

Certified as sterile with a sterility assurance level (SAL) of 10^{-6}