



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 12x75 mm Culture Tubes, Polypropylene, 125 per Bag, Sterile

Fluid Measure	mL	Freezable	Yes
Fluid Volume	5	Freeze Temp	-80°C
Height (mm)	75	Autoclavable	121°C 15min
Width (mm)	12	Sterile	Yes
Length (mm)		Sterile Processing	Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶
Pack Qty	125	Centrifugation	Yes
Case Qty	1000	Centrifugation Speed (RCF)	3000
Package Type	Bag	Shelf Life (Months)	60
Material	Polypropylene	Certification	ISO 9001
Cap Material		Regulatory	U.S. FDA 21 CFR

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing if applicable:

ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)



This Technical Data Sheet is provided as a service and for documentation only. Labcon North America retains all rights to the information, procedures, and results of the material contained hereon as proprietary information. Any copying and distribution of this technical data without prior written authorization is prohibited by law. The use of this information and documentation does not discharge the user from their responsibility to test and use this product appropriately. This item is for use solely by laboratory professionals. This document is automatically generated and is © Copyright 05-Nov-2025 . All Rights Reserved. Technical Data Sheet Version 1 Rev 1/2023

DOMINIQUE DUTSCHER S.P.A.