



Dimensions are inches. Do not scale



Image is not actual size. Colors may differ.

SuperClear® 17x100 mm Culture Tubes with Dual Position Caps, Polypropylene, 125 per Bag, Sterile

| | | | |
|---------------|---------------|----------------------------|---|
| Fluid Measure | mL | Freezable | Yes |
| Fluid Volume | 14 | Freeze Temp | -80°C |
| Height (mm) | 100 | Autoclavable | 121°C 15min |
| Width (mm) | 17 | Sterile | Yes |
| Length (mm) | | Sterile Processing | Certified as sterile with a sterility assurance level (SAL) of 10 ⁻⁶ |
| Pack Qty | 125 | Centrifugation | Yes |
| Case Qty | 1000 | Centrifugation Speed (RCF) | 3000 |
| Package Type | Bag | Shelf Life (Months) | 36 |
| Material | Polypropylene | Certification | |
| Cap Material | Polyethylene | Regulatory | |

Each manufacturing lot is sampled and tested in accordance with our Standard Operating Procedures.

Additional Testing: ATP, RNase/DNase and Endotoxin Free (non-pyrogenic)

