

## Description

- Blood cell counters

## Part number

- 952507

Pieces by sales : 1



---

## Description

- Electronic counters and memory devices for microscopic blood testing
- The electronic memory devices specify blood counting.
- The microscope image can be observed permanently, the numbers can be typed in easily. Every device comes with an instruction manual.
- Advantages :
  - Power supply operation
  - Rechargeable battery operation
  - Software (download) with computer analysis
  - Big TFT color display
  - Storage of last counting on device
- Interfaces : Data transfer – USB mini
- Power supply
  - Power supply operation – USB adapter 5VDC / 300mA
  - Rechargeable – USB power bank external
  - Battery operation
- Functional keys : Labelling On/Off, -/Eval, +, Go On, Corr. Clear
- Counting keys
  - Labelling : 40345010 – Eos, Baso, Promy, Myelo, Div, Juve, Stab, Seg, Lym, Mono
  - 40345015 – At Ly, My Bl, Met My, Myelo, Normo, Lym, Seg, Mono, Plas, Tox Gr, Stab, Eos, Baso, Div/Reti, Div/Ery
- Programs
  - Counting steps – 10-1000 in 10 steps
  - Programs – LEUKO, MYELO/RBC, RETI, FREE

## Technical data

Dimensions : 130 x 235 x 35 mm

Number of functional keys : 5

Colour : blue

Counting limit : max. 1000

Display : TFT color display 3,5'