

UNITED POLY ENGINEERING PVT. LTD.

D-13/3, OKHLA INDUSTRIAL AREA

PHASE – II, NEW DELHI, INDIA

www.unitedpoly.com



SURGICAL TRAY WITHOUT COVER, STAINLESS STEEL (AISI 304)



AUSTENITIC STAINLESS STEEL

CODE
UPL-025

DIM (MM)
355x254x50

Technical Specification:

Basic item of equipment for carrying surgical instruments, dressing materials, injection materials, etc.

Can be sterilized in a steam sterilizer (autoclave) for sterile purposes.

Material: Austenitic Stainless Steel (AISI 304) with a smooth surface.

Austenitic Stainless Steel composition: approximately 8 to 10% nickel, 18 to 20% chromium.

Thickness: 0.6 mm.

Size tolerance: approximately ± 3 mm.