

UNITED POLY ENGINEERING PVT. LTD.

D-13/3, OKHLA INDUSTRIAL AREA
PHASE – II, NEW DELHI, INDIA

www.unitedpoly.com



INSTRUMENT TRAY, STAINLESS STEEL



STAINLESS STEEL AISI 304 GRADE

CODE
UPL -033

DIM (MM)
313 x 218 x 31

TECHNICAL SPECIFICATION:

Product: Stainless Steel AISI 304 Grade Instrument Tray

Description: The stainless steel AISI 304 grade instrument tray (UPL-033) is designed as basic equipment for nursing and surgical care. It is suitable for use in medical environments and can be sterilized in a steam sterilizer or hot air sterilizer for sterile purposes.

Material: The instrument tray is constructed from stainless steel with a smooth surface. The stainless steel composition consists of approximately 8% nickel and 18% chromium, providing durability and corrosion resistance.

Specifications:

- Thickness: 0.6 mm
- Finish: Available in bright or satin finish options
- Tolerance: Dimensions are held to a tolerance of ± 2 mm